

## Optical Halogen Lamp 1003188

### Instruction Sheet

10/13 ALF



- 1 Hole diaphragm
- 2 Rod
- 3 Condenser
- 4 Retaining ring
- 5 Lamp housing
- 6 Air slits

### 1. Safety instructions

Caution. Lamps get hot when they are switched on for a long period of time.

- When handling the lamp during the experiment, use a cloth or other protection.
- After the experiment leave the lamp to cool.
- Do not cover the air slits.
- Do not supply the lamp with an operating voltage in excess of 12 V.

### 2. Description

The optical halogen lamp is an ultra-bright light source for experiments on optical bench and for projection.

It consists of a metal housing with a built-in condensing lens and integrated cooling fan. The knob on the rear for moving the lamp bulb makes it possible to obtain a parallel light beam. To reduce

the intensity of the light beam, a hole diaphragm can be attached to the condensing lens by means of a retaining ring. The handling rod can be unscrewed if necessary.

### 3. Accessories

Spare halogen lamp, 12 V, 50 W	1002837
--------------------------------	---------

### 4. Technical data

Halogen lamp:	12 V, 50 W
Connection:	4-mm safety jacks
Condenser focal length:	75 mm
Condenser diameter:	45 mm
Rod:	120 mm x 10 mm diam.
Dimensions:	approx. 190x125x110 mm <sup>3</sup>

## 5. Operation

### 5.1 General notes

Suitable voltage sources for operation of the halogen lamp are, for example

Transformer 12 V, 60 W (230 V) 1000593  
or

Transformer 12 V, 60 W (115 V) 1006780

- Do not allow the lamp to suffer mechanical shocks.

### 5.2 Producing a parallel light beam

- Set up the optical lamp on the optical bench or clamp it in a barrel foot and connect it to the voltage source.
- Set up the screen at the required distance from the lamp.
- Loosen the locking screw of the bulb movement knob on the rear of the lamp.



Fig. 1 Rear side of the lamp: 1 fan, 2 locking screw, 3 bulb movement knob, 4 sockets

- Slide the bulb movement knob forward or backward until a sharp image of the bulb filament is obtained on the screen.
- Re-tighten the locking screw.
- To reduce the intensity of the light beam, attach the hole diaphragm to the condenser lens using the retaining ring.

### 5.3 Changing the bulb

Any deposits of fat from the skin on a halogen lamp bulb cause the glass to become fogged and significantly reduce the life of the bulb.

- Do not touch the glass bulb of the halogen lamp with fingers.

- Remove the four securing screws and remove the lid of the casing.
- Take out the defective halogen bulb.
- Holding the 12 V, 50 W replacement bulb with a cloth or other protective material, push the connecting wires into the socket.
- Screw the lid back on.




Fig. 2 Changing the bulb: 1 securing screws, 2 lid of the casing, 3 socket

### 5.4 Setting up on the optical bench

- Screw the handling rod into the base of the optical lamp.
- Fix the optical lamp into the rider on the optical bench.



Fig. 3 Setting up on the optical bench