



  
C-Celia®

Emergency Obstetrics  
Simulation Suite

Fetal Extraction Simulator

Product Number: OB-500-2000

1022559

## Table of Contents

<b>SECTION 1: IMPORTANT INFORMATION</b> .....	3
<b>SECTION 2: Introduction</b> .....	4
Skills Taught .....	4
Items Included .....	4
<b>SECTION 3: Unpacking and Assembly</b> .....	5
<b>SECTION 4: Directions for Use</b> .....	5
Preparing the Simulator .....	5
Loading the Baby .....	6
<b>SECTION 5: Repairing the Simulator After Use</b> .....	7
<b>SECTION 6: Care and Maintenance</b> .....	7
Replacement Parts .....	7
Customer Support .....	8

# Important Information

## GENERAL WARNINGS

---

**Warning: A poorly lubricated baby and placenta may seriously damage the simulator.**

Only use the lubricant supplied with this product. Use of non-recommended fluid products may weaken or damage the simulator.

Lubrication is key to the usage of this product; inadequate lubrication could lead to tears or other damage to silicone parts. Apply at least 2-3 squirts of lubrication to the following areas:

- Inside the abdominal cavity
  - The baby's head, shoulders, body, limbs and placenta.
- Ensure the baby is well covered.

It is vitally important that after each training session all lubrication is cleaned with a warm damp cloth. You can also use baby wipes to wipe away any excess lubricant.

**Warning: Sharp objects, excessive force, or using non-OEI products can scratch, tear or damage the simulator.**

Gloves should always be worn, and all jewelry items should be removed. Sharp surfaces such as finger nails and rings can sometimes lead to damage of silicone parts.

## WARRANTY INFORMATION

---

WARRANTIES. EXCEPT AS OTHERWISE SET FORTH HEREIN, OPERATIVE EXPERIENCE, INC. ("OPERATIVE EXPERIENCE") PRODUCTS ARE SOLD ONLY WITH THE 1-YEAR LIMITED WARRANTY DESCRIBED BELOW. THIRD-PARTY PRODUCTS ARE EXPRESSLY EXCLUDED FROM THE BELOW LIMITED WARRANTY, AND ARE SUBJECT ONLY TO THOSE WARRANTIES EXTENDED BY THE APPLICABLE THIRD PARTY MANUFACTURERS OF SUCH PRODUCTS, IF ANY. EXCEPT WHERE PROHIBITED BY LAW OR AS OTHERWISE PROVIDED HEREIN, OPERATIVE EXPERIENCE HEREBY EXPRESSLY DISCLAIMS ALL WARRANTIES EITHER EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, NON-INFRINGEMENT, AND ANY WARRANTIES ARISING FROM A COURSE OF DEALING, COURSE OF PERFORMANCE, OR USAGE OF TRADE. OPERATIVE EXPERIENCE PERSONNEL, SUPPLIERS, AND LICENSEES ARE NOT AUTHORIZED TO ALTER THESE WARRANTIES.

Operative Experience warrants that its products ("Products") are free from defects in materials and workmanship, subject to specific limitations and conditions set forth below, for a period of one (1) year from the date of the original purchase by the first end-user (the "Warranty Period"). Claims of defect brought after the Warranty Period will not be covered by this limited warranty. In order to bring a claim under this limited warranty, the buyer must notify Operative Experience by email to [info@operativeexperience.com](mailto:info@operativeexperience.com), describing the defect in reasonable detail. If Operative Experience determines that such defect may be covered by this limited warranty, Operative Experience will issue a return material authorization ("RMA") to the buyer, and the buyer shall return the Product in question to Operative Experience using the instructions provided by Operative Experience. If Operative Experience, in its reasonable discretion, determines that a Product has a defect in materials or workmanship, Operative Experience shall, at its option, repair or replace the defective Product at no cost to the buyer. If identical materials are not available at the time of repair or replacement, Operative Experience may, at its option, substitute materials of comparable quality or refund buyers money. The Products are intended only for use in medical training, and this limited warranty shall not cover any Product used for commercial or other non-medical training purposes. Any alteration by buyer of a Product shall render this limited warranty void with respect to the altered Product. This limited warranty extends only to the original buyer and shall not cover defects or damages resulting from (i) accident, misuse, abuse, tampering, neglect, unusual physical stress, alteration, or improper operation or handling of the Product; (ii) incorrect or inadequate maintenance, cleaning, or care of the Product; or (iii) normal wear and tear of the Product. Operative Experience will not be responsible under this limited warranty or otherwise for defects or damage caused by any ancillary equipment not furnished by Operative Experience which is attached to or used in connection with the Products, and all such equipment is expressly excluded from this limited warranty. This limited warranty does not extend to consumable parts, such as batteries, or any third-party products that buyer may use or purchase in conjunction with the Products.

# INTRODUCTION

The Fetal Extraction Simulator is designed for clinicians and learners to experience and practice difficult deliveries such as shoulder dystocia, breech, nuchal cord, prolapse arm and use of suction assist and forceps.

Open (post c-section) abdomen provides a true-to-life difficult delivery operative experience

Low-transverse incision and uterus with a low transverse Hysterotomy enable various extraction maneuvers and techniques including suction assist and forceps

## Skills Taught

- Normal and high-risk deliveries
- Instrument assisted delivery
- Delivery of the placenta
- Teamwork and human factors training in emergency Caesarean section

## Items Included

<b>1</b> Simulator Base Assembly Part Number: OB-500-2000	Qty: 1	<b>4</b> Lubricant, 16oz/473ml bottle Part Number: 220007	Qty: 1
<b>2</b> Small Baby Part Number: OB-050-2925	Qty: 2	<b>5</b> Spare Hardware Kit Part Number: OB-300-2400	Qty: 1
<b>3</b> Silicone Glue Part Number: OB-050-0998	Qty: 1		



# Unpacking and Assembly

---

**NOTE: Use caution when opening the shipping carton, as this is preferred storage container. Damage to the container could result in an inability to store it properly, leading to damage of the simulator and accessories.**

C-Celia comes in a specially designed carton with cut-outs in the internal foam that provide for placement of all parts and accessories. Keep all original shipping materials, including boxes.

When unpacking C-Celia for the first time, careful use of a box cutter protects both the packaging and the product.

It is recommended that you reserve a space approximately 2 feet (24 inches) by 6 feet (72 inches) into which you can place the contents of shipping carton. Remove each item in the carton and verify, on the shipping paperwork, that it is the correct part.

# Directions for Use

---

## Preparing the Simulator



Prior to simulating a delivery, the simulator needs to be fully lubricated before use.



Using the lube provided, generously apply lube to the abdominal cavity using your hands.

# Loading the Baby

To load the baby and placenta:



1 Generously lubricate the baby and placenta.



2 Open the back plate



3 Remove inner door



4 Insert the baby and placenta into the desired position



5 Replace inner door



6 Once the baby is positioned, close and secure the back plate

# Repairing the Simulator After Use

---

## Repairing the Umbilical Cord

### Materials needed

- Paper towels
- Popsicle Sticks
- Acetone
- Silicone Glue
- Pins

### Gluing the cut back together

1. Using paper towels clean off all lubricant from the umbilical cord.
2. Clean the entire area on the umbilical cord that needs to repair with Acetone. Let dry
3. Using a popsicle stick or a gloved hand apply a thin layer of Silicone Glue to the two cut ends.

4. Using the pins put 2-3 through both sides of the cord to hold in place and let dry.
5. Allow the Silicone Glue to dry for about 25 minutes or 5 mins with a hairdryer or heat gun.
6. Remove the pins

# Care and Maintenance

---

### Surface

Use only on clean surface. Avoid felt tipped markers, ink pens, acetone, iodine or other staining products and avoid placing the module/manikin on newsprint or inked lines of any kind.

### Lubrication

Cottonseed Oil should be used liberally to lubricate both the baby and placenta and the inside of the abdominal cavity.

### Storage

All fluid should be removed from the abdominal cavity before storage. All surface areas should be cleaned and dried before storage.

### Cleaning

Simulators should be wiped clean with a solution of mild soap and water. Clean with mild soap and water; do not submerge simulator, module or parts in cleaning fluids or water. The entire model can be rinsed using a hose. Parts should be drained and air-dried thoroughly before storage and disinfected when needed.

## Replacement Parts

Description	Part Number
Shell Door	OB-100-2005
Lubricant	220007
OB Spare Parts Kit	OB-300-2400



MORE INFORMATION AVAILABLE ONLINE AT [3BSCIENTIFIC.COM](http://3BSCIENTIFIC.COM)



**3B Scientific GmbH** • Ludwig-Erhard-Straße 20 • 20459 Hamburg • Germany • [3bscientific.com](http://3bscientific.com)  
Phone: +49 40 73966-0 • Fax: +49 40 73966-100 • E-mail: [info@3bscientific.com](mailto:info@3bscientific.com)