
Transformator Kern D 1022664 Polschuhe D 1000978

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10/20 JH/ GH

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Instruction sheet

10/20 JH/ GH

- 1 Transformer core (U-core)
- 2 Yoke (I-core)
- 3 Clamps for clamping the yoke or pole pieces firmly to the core.
- 4 Pair of drilled pole pieces

1. Safety instructions

- Polished surfaces should be kept free of dirt and grease.
- Transformer core, yoke and pole pieces should not be exposed to moisture.
- For transporting the U-core and yoke (I-core), make sure the clamps are firmly secured.
- During the experiment, the yoke or pole pieces should be firmly secured using the clamps.

The transformer core with yoke is a tape wound core with two holes drilled for securing pole pieces or yoke with the aid of clamps. For better stability the U-core is mounted on a plastic base.

Cross section of core:	38 mm x 40 mm
Height including yoke:	170 mm approx.
Width:	150 mm approx.
Material:	Iron, tape wound core
Weight:	5.3 kg approx.

2. Description, technical data

2.1 Transformer core

The transformer core with yoke and clamps are provided for use in conjunction with the accessories listed under item 4 of the instructions for assembling the demountable transformer.

- 2.1.1 Scope of delivery
- 1 Transformer core
 - 1 Yoke
 - 2 Clamps (pair)

2.2 Pair of drilled pole pieces

The pole pieces are required for experiments on electromagnetism where a well-defined air gap is necessary (e.g. Waltenhofen's pendulum or investigations of paramagnetic and diamagnetic samples).

Pole pieces made of soft iron each with one plane and one conical end. The pole pieces have a bore drilled through the middle of them.

Cross section of core:	40 mm x 40 mm
Bore at conical end:	5 mm
Bore at plane end:	12 mm
Material:	Soft iron
Weight:	1.7 kg approx.

3. Operation

- The safety instructions for the coils must be observed at all times.
- Mount the primary and secondary coils on the core.
- Lay the polished side of the yoke or pole pieces on top of the U-core.
- Attach the clamps.
- Firmly secure the yoke or pole pieces with the clamps.

4. Accessories for demountable transformer

Mains coil D @230V 220V	1000987
Mains coil D @115V	1000986
Low Voltage Coil D	1000985
Coil D (600 turns)	1000988
Coil D (1200 turns)	1000989
Coil D (6000 turns)	1000990
Fusion ring	1000980
High current coil with 5 turns for spot welding	1000981
High current coil D for nail fusing experiment	1000984
High voltage coil D including 2 horn shaped electrodes	1000991
Waltenhofen's pendulum	1000993

5. Cleaning, maintenance and storage

- Polished surfaces should be kept free of dirt and grease.
- Store in dry conditions.
- Remove any rust with fine steel wool or sandpaper.

6. Example experiments

Fig.1 Demountable transformer

Fig. 2 Waltenhofen's pendulum



Fig. 3 Spark discharge along hornshaped electrodes

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